Made in Belgium

### REF

Wheeleo-v2	Wheeleo White
Wheeleo-v2/9005	Wheeleo Black
Wheeleo-v2/3003	Wheeleo Red
Wheeleo-v2/5009	Wheeleo Blue
Wheeleo-v2/5021	Wheeleo Green

### Product Name

NOTICE D'UTILISATION	FR	2
GEBRUIKAANWIJZING	NL	7
GEBRAUCHSANWEISUNG	DE	12
USER MANUAL	EN	17
MANUALE D'USO	IT	22
MANUAL DE INSTRUCCIONES	ES	27
MANUAL DE UTILIZAÇÃO	PT	32
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	GR	37
تعليمات الاستخدام	AR	42

# Wheeleo®

One hand for autonomy

## 1. Généralités

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**#onehandforautonomy**

## 1. General

InnoRehab would like to thank you for your purchase of the Wheeleo<sup>®</sup>. This user manual contains information on the correct use of the product and important safety advice.

Please read the information contained in this manual carefully before using the walking cane. Please consult a healthcare professional in order to ensure that the cane is suitable for your needs.

Please contact your stockist for additional information, advice or assistance.

## 2. Intended standard use

The Wheeleo<sup>®</sup> is a hand-held device to aid mobility for people with walking difficulties.

## 3. Contraindications

Not to be used by people whose weight exceeds 120 kg.  
Not to be used for any purpose other than that specified in these instructions.

EN

## 4. Declaration of conformity

Wheeleo<sup>®</sup> is a walking aid medical device that complies with the applicable provisions described in the Medical Equipment Regulation (EU) 2017/745.

## 5. Warnings

Read the whole of this manual carefully before using the Wheeleo<sup>®</sup>.

Check that the Wheeleo<sup>®</sup> is supplied complete. If it is to be stored for a lengthy period, keep the Wheeleo<sup>®</sup> wrapped and away from heat and dampness.

The service life of the Wheeleo<sup>®</sup> is determined by the wear on its component parts.

The Wheeleo is designed for use on flat ground and relatively smooth ground. The castors must be able to roll freely.

InnoRehab does not accept any liability for the consequences of improper use of the product or for any modifications made to the Wheeleo<sup>®</sup> structure.


In the event of a change in performance or malfunction of the Wheeleo, InnoRehab should be contacted and the product should not be used again until the problem is resolved.

## 6. Description

### a. Wheeleo<sup>®</sup>

The Wheeleo<sup>®</sup> is made of steel with a paint finish. It has 4 100mm castors enabling it to turn through 360°. It incorporates a mechanism for rapid adjustment of the handle height. The

height of the cane can be adjusted in 2 cm steps between 77,5 cm and 101,5 cm. The Wheeleo® allows the user to walk quickly, stably and comfortably. The maximum weight of the user is 120 kg.

N°	Description
1	Handle
2	Sliding tube
3	Anneau finition
4	Structure of the foot
5	Index Plunger
6	Castor

**b. Adjustment of the height**

To adjust the height of the handle, grip the stopper in three fingers and pull out, then use the other hand to slide the sliding tube over the central tube to connect with a higher or lower adjustment hole. Release the stopper in order to allow it to re-engage. Ensure that the holes in the sliding tube are in line with the adjustment hole. If you are not able to do this alone, ask someone to help you. A healthcare professional can advise on the ideal height for your requirements.

For a more occasional balance aid, an adjustment higher than the hip is recommended. For a 'suportive' type of aid, adjusting the handle to the same height as the bend in your wrist may be more comfortable. In all cases, a professional can advise you on the ideal height for your needs.

Setting


Pull in the axis of the adjustment plug !

Do not unscrew !

**c. Use in the right or left hand**

No adjustment to the handle position is necessary in order to use the Wheeleo<sup>®</sup> in the right or left hand. Just turn the Wheeleo<sup>®</sup> round. The handle is symmetrical. On the other hand, the Wheeleo<sup>®</sup> is not completely symmetrical. It has a narrow side and a wider side. The wider side and the indexing button must be closest to the user. The green dot must be visible to the user.

EN

Left hand

Right hand

## 7. Before use

- Check that the structure of the Wheeleo® is not damaged in order to ensure it is completely safe to use.
- Check that the castors are firmly fixed and roll smoothly. The castors must also swivel freely through 360°.
- It is advisable to be accompanied by a professional when you first use the Wheeleo® in order to maximise the benefit.
- It is advisable to modify the location where the Wheeleo® is used in order to reduce the risk of falls. The floor must be as smooth as possible. It is advisable to remove rugs.
- Check that it is the correct height and that the adjustment mechanism is correctly positioned.

## 8. Warnings for correct use

- Do not use if your weight exceeds 120 kg.
- Risk of tipping: keep the Wheeleo® facing in the right direction (section 6.c) and leave the casters free to move. The Wheeleo is not suitable for floors with excessively uneven surfaces that prevent the castors from rolling smoothly. For example (non-exhaustive list) carpets, gravel, paving stones, electric cables, grass, etc.
- Danger of tipping over when all 4 wheels of the Wheeleo® are off the ground. Do not tilt the Wheeleo® more than 5° outwards or 2° inwards.
- Avoid using the Wheeleo® on wet surfaces. Water can damage the bearings and prevent the wheels from turning.
- Avoid using the Wheeleo® outdoors, on dirty surfaces (sand, dust, etc.). Dirt can also prevent the wheels from turning properly.
- If use of the Wheeleo® causes problems (allergic reaction or pain), stop using it and contact a healthcare professional.
- Risk of trapping the fingers: do not insert fingers into the holes in the handle.
- Any serious incident with the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

## 9. Maintenance

The Wheeleo<sup>®</sup> was verified and certified with the label “CE – Medical Device Class 1” when launched on the market. For the user’s safety, it is recommended that you have the device checked by the manufacturer or an approved laboratory every 2 years. In the event of repair, only use original equipment spare parts or accessories.

## 10. Cleaning and disinfection

- a. **Cleaning:** Use a cloth to remove dust. If necessary, add a small amount of detergent to remove stubborn dirt. Dry it completely before use.
- b. **Disinfection:** If you need to disinfect the equipment, use a disinfectant cleaner such as hydroalcoholic gel.

**Warning: Do not use acids, alkalis or solvents such as acetone.**

## 11. Disposal

Do not dispose of the device in normal refuse bins. Dispose of the Wheeleo<sup>®</sup> at an appropriate recycling centre.

EN

## Product Registration :

Online: <http://www.wheeleo.be/registration>

Or by sending the document on page 47 completed

**Connect to the Wheeleo<sup>®</sup> community to share your experiences and get advice from other users!**

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A registered Trademark of InnoRehab SRL