Wheeleo®

One hand for autonomy

The first One-handed Rollator

REF

Wheeleo-v2

Wheeleo-v2/5009

Wheeleo-v2/9005

Product Name

Wheeleo White

Wheeleo Blue

Wheeleo Black

Produced with love in Belgium

EU PATENT: EP4380532

More Info on www.wheeleo.eu/en/patents

NOTICE D'UTILISATION	FR	2
GEBRUIKAANWIJZING	NL	10
GEBRAUCHSANWEISUNG	DE	15
USER MANUAL	EN	20
MANUALE D'USO	IT	25
MANUAL DE INSTRUCCIONES	ES	30
MANUAL DE UTILIZAÇÃO	PT	35
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	GR	40
تعليمات الاستخدام	AR	45
BRUGSANVISNING	DK	50

More languages on

<https://www.wheeleo.eu/operating-instructions/>

Senior and modern... Incompatible ?

Wheeleo®
One hand for your autonomy

Wheeleo®,
Your new-generation
walking companion that
ensures stability and
safety throughout
your day.

Clinical
trial

Test me !

WHAT DO USERS THINK?

"It's the **best indoor rollator**. The Wheeleo® is compact and space-saving compared to my classic rollator. It moves **effortlessly** in all directions."

Marie-Paule C.
93 year old user

"I **feel safer** with the Wheeleo®. If I didn't have the Wheeleo®, I would lose my autonomy, I would surely be in my chair not knowing how to move around."

Vincent Durant
37 year old hemiplegic

Smart design

Compact. Convertible. Elegant.

Our design allows **easy transportation** and has a **modern look**. No stigma, just a next-generation mobility aid

One hand

Experience full support, without compromising your freedom. Always one hand free for daily activities.

4 points of contact

4 wheels that keep you **connected to the ground at all times**.

Experience the stability of a rollator's 4 legs, in the shape of a cane.

Dynamic braking

The Wheelleo® has no brakes, and more importantly, doesn't need brakes. The device is **easily controlled**, enabling the person to stop and walk at a desired pace.

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#onehandforautonomy

1. General

InnoRehab would like to thank you for your purchase of the Whee leo®. This user manual contains information on the correct use of the product and important safety advice.

Please read the information contained in this manual carefully before using the walking cane. Please consult a healthcare professional in order to ensure that the cane is suitable for your needs.

Please contact your stockist for additional information, advice or assistance.

2. Intended standard use

The Whee leo® is a hand-held device to aid mobility for people with walking difficulties.

3. Contraindications

Not to be used by people whose weight exceeds 120 kg.
Not to be used for any purpose other than that specified in these instructions.

4. Declaration of conformity

Whee leo® is a walking aid medical device that complies with the applicable provisions described in the Medical Equipment Regulation (EU) 2017/745.

5. Warnings

Read the whole of this manual carefully before using the Whee leo®.

Check that the Whee leo® is supplied complete. If it is to be stored for a lengthy period, keep the Whee leo® wrapped and away from heat and dampness.

The service life of the Whee leo® is determined by the wear on its component parts.

The Whee leo is designed for use on flat ground and relatively smooth ground. The castors must be able to roll freely.

InnoRehab does not accept any liability for the consequences of improper use of the product, a lack of maintenance, or for any modifications made to the Whee leo® structure.


In the event of a change in performance or malfunction of the Whee leo, InnoRehab should be contacted and the product should not be used again until the problem is resolved.

6. Description

a. Whee leo®

The Whee leo® is made of steel with a paint finish. It has 4 100mm castors enabling it to turn though 360°. It incorporates a mechanism for rapid adjustment of the handle height. The

height of the cane can be adjusted in 2 cm steps between 77,5 cm and 101,5 cm. The Wheeleo® allows the user to walk quickly, stably and comfortably. The maximum weight of the user is 120 kg.

N°	Description
1	Handle
2	Sliding tube
3	Anneau finition
4	Structure of the foot
5	Index Plunger
6	Castor


EN

b. Adjustment of the height

To adjust the height of the handle, grip the stopper in three fingers and pull out, then use the other hand to slide the sliding tube over the central tube to connect with a higher or lower adjustment hole. Release the stopper in order to allow it to re-engage. Ensure that the holes in the sliding tube are in line with the adjustment hole. If you are not able to do this alone, ask someone to help you. A healthcare professional can advise on the ideal height for your requirements.

For a more occasional balance aid, an adjustment higher than the hip is recommended. For a 'suportive' type of aid, adjusting the handle to the same height as the bend in your wrist may be more comfortable. In all cases, a professional can advise you on the ideal height for your needs.

Setting


Pull in the axis of the adjustment plug !

Do not unscrew !

c. Use in the right or left hand

No adjustment to the handle position is necessary in order to use the Wheeleo® in the right or left hand. Just turn the Wheeleo® round. The handle is symmetrical. On the other hand, the Wheeleo® is not completely symmetrical. It has a narrow side and a wider side. The wider side and the indexing button must be closest to the user. The green dot must be visible to the user.

Left hand

Right hand

7. Before use

- Check that the structure of the Wheeleo® is not damaged in order to ensure it is completely safe to use.
- Check that the castors are firmly fixed and roll smoothly. The castors must also swivel freely through 360°.
- It is advisable to be accompanied by a professional when you first use the Wheeleo® in order to maximise the benefit.
- It is advisable to modify the location where the Wheeleo® is used in order to reduce the risk of falls. The floor must be as smooth as possible. It is advisable to remove rugs.
- Check that it is the correct height and that the adjustment mechanism is correctly positioned.

8. Warnings for correct use

- Do not use if your weight exceeds 120 kg.
- Risk of tipping: keep the Wheeleo® facing in the right direction (section 6.c) and leave the casters free to move. The Wheeleo is not suitable for floors with excessively uneven surfaces that prevent the castors from rolling smoothly. For example (non-exhaustive list) carpets, gravel, paving stones, electric cables, grass, etc.
- Danger of tipping over when all 4 wheels of the Wheeleo® are off the ground. Do not tilt the Wheeleo® more than 5° outwards or 2° inwards.
- Avoid using the Wheeleo® on wet surfaces. Water can damage the bearings and prevent the wheels from turning.
- Avoid using the Wheeleo® on dirty surfaces (sand, dust, etc.). Dirt can also prevent the wheels from turning properly.
- If use of the Wheeleo® causes problems (allergic reaction or pain), stop using it and contact a healthcare professional.
- Risk of trapping the fingers: do not insert fingers into the holes in the handle.
- Any serious incident with the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

9. Maintenance

The Wheeleo® was verified and certified with the label “CE – Medical Device Class 1” when launched on the market. For the user’s safety, it is recommended that you have the device checked by the manufacturer or an approved laboratory every 2 years. In the event of repair, only use original equipment spare parts or accessories.

On a regular basis and especially if you use the Wheeleo® on uneven ground, which can cause vibrations, check that the wheels are still screwed tight. If necessary, tighten them again (14 mm wrench). You can add threadlocker to reinforce the tightening. For more information on maintenance, visit <http://www.wheeleo.be/registration>

10. Cleaning and disinfection

- a. **Cleaning:** Use a cloth to remove dust. If necessary, add a small amount of detergent to remove stubborn dirt. Dry it completely before use.
- b. **Disinfection:** If you need to disinfect the equipment, use a disinfectant cleaner such as hydroalcoholic gel.

Warning: Do not use acids, alkalis or solvents such as acetone.

11. Logo

Manufacturer's name and address

Year of manufacture of the device

Device reference number

Device serial number

Product Registration and Maintenance :

Online: <https://www.wheeleo.eu/registration>

For registration, you can send the document on page 55 completed

Connect to the Wheeleo® community to share your experiences and get advice from other users!

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A registered Trademark by InnoRehab SRL